




VERSYS® POROUS REVISION STEMS

Multiple Options
In Your Hands

VerSys®
HIP SYSTEM

STANDARD

REVISION SYSTEMS







STANDARD BEADED FULLCOAT REVISION STEMS

LM Body • CoCrMo • Revision Offset

VerSys[®]
HIP SYSTEM




	Item Number		Offset at +0mm Head Ctr (mm)	Diameter (mm)	Stem Length	Anteversion	Strength Enhancement Feature
	Bowed Left	Bowed Right					
 8" Bowed Revision Stems	00-7843-013-81	00-7843-013-82	42	13.5	200mm	8° Head/Neck 8° Proximal Body	Sizes 13.5 & 15 only
	00-7843-015-81	00-7843-015-82	45	15			
	00-7843-016-81	00-7843-016-82	45	16.5			
	00-7843-018-81	00-7843-018-82	48	18			
	00-7843-019-81	00-7843-019-82	48	19.5			
	00-7843-021-81	00-7843-021-82	48	21			
	00-7843-022-81	00-7843-022-82	48	22.5			
 10" Bowed Revision Stems	00-7843-013-11	00-7843-013-12	42	13.5	250mm	8° Head/Neck 8° Proximal Body	Sizes 13.5 & 15 only
	00-7843-015-11	00-7843-015-12	45	15			
	00-7843-016-11	00-7843-016-12	45	16.5			
	00-7843-018-11	00-7843-018-12	48	18			
	00-7843-019-11	00-7843-019-12	48	19.5			
	00-7843-021-11	00-7843-021-12	48	21			
	00-7843-022-11	00-7843-022-12	48	22.5			
 7.5" Straight Revision Stems	00-7843-012-18		39	12	190mm	No Anteversion	Sizes 12, 13.5, & 15 only
	00-7843-013-18		42	13.5			
	00-7843-015-18		45	15			
	00-7843-016-18		45	16.5			
	00-7843-018-18		48	18			
	00-7843-019-18		48	19.5			
	00-7843-021-18		48	21			
00-7843-022-18		48	22.5				
 10" Straight Revision Stems	00-7843-015-24		45	15	250mm	No Anteversion	Size 15 only
	00-7843-016-24		45	16.5			
	00-7843-018-24		48	18			
	00-7843-019-24		48	19.5			
	00-7843-021-24		48	21			
	00-7843-022-24		48	22.5			

BEADED FULLCOAT CALCAR REVISION STEMS


LM Body • CoCrMo • Revision Offset • +10mm and +20mm Calcar Build-ups on Stem

VerSys[®]
HIP SYSTEM



	Item Number		Offset at +0mm Head Ctr (mm)	Diameter (mm)	Stem Length	Anteversion	Strength Enhancement Feature
	Bowed Left	Bowed Right					
10" Bowed Revision +10 Calcar Stems	00-7843-013-91	00-7843-013-92	42	13.5	240mm	8° Head/Neck 8° Proximal Body	Sizes 13.5 & 15 only
	00-7843-015-91	00-7843-015-92	45	15			
	00-7843-016-91	00-7843-016-92	45	16.5			
	00-7843-018-91	00-7843-018-92	48	18			
	00-7843-019-91	00-7843-019-92	48	19.5			
	00-7843-021-91	00-7843-021-92	48	21			
	00-7843-022-91	00-7843-022-92	48	22.5			
10" Bowed Revision +20 Calcar Stems	00-7843-015-71	00-7843-015-72	45	15	230mm	8° Head/Neck 8° Proximal Body	Size 15 only
	00-7843-016-71	00-7843-016-72	45	16.5			
	00-7843-018-71	00-7843-018-72	48	18			
	00-7843-019-71	00-7843-019-72	48	19.5			
	00-7843-021-71	00-7843-021-72	48	21			
	00-7843-022-71	00-7843-022-72	48	22.5			
7.5" Straight Revision +20 Calcar Stems	00-7843-015-17		45	15	167mm	No Anteversion	Size 15 only
	00-7843-016-17		45	16.5			
	00-7843-018-17		48	18			
	00-7843-019-17		48	19.5			
	00-7843-021-17		48	21			
	00-7843-022-17		48	22.5			
7.5" Straight Revision +10 Calcar Stems	00-7843-013-07		42	13.5	177mm	No Anteversion	Sizes 13.5 & 15 only
	00-7843-015-07		45	15			
	00-7843-016-07		45	16.5			
	00-7843-018-07		48	18			
	00-7843-019-07		48	19.5			
	00-7843-021-07		48	21			
	00-7843-022-07		48	22.5			

CALCAR





zimmer
Confidence in your hands™